


# INSTRUCTION MANUAL LINOMAT 5





### Table of content

1	Introduction	2
1.1	CAMAG Linomat 5	2
1.2	Precautions	2
1.3	Parts supplied	4
1.4	Spare parts	4
2	Installation	5
2.1	Instrument	5
2.2	Installation environment	6
2.3	Conditions for the installation	7
2.4	Gas supply	7
2.5	Mains connection	7
2.6	Connection to the PC	8
3	Getting started	8
3.1	Automatic instrument monitoring	8
3.2	Checking system functions	8
3.3	Inserting a dosage syringe	9
3.4	Setting the application position	10
3.5	Loading a TLC plate or foil	10
3.6	Waste plate	11
4	Manual operation	12
4.1	Display and function of keys	12
4.2	Stand-alone operation	13
4.3	Aborting an application in progress	14
4.4	Resetting the instrument	14
4.5	Docking the syringe	14
5	Maintenance and trouble shooting	15
5.1	General maintenance	15
5.2	The Setup dialog	15
Set ba	aud rate	15
Set dis	splay language	15
Set th	e LCD contrast	16
5.3	Handling of the HAMILTON-syringe	16
5.4	Adjusting the spray gas	16
5.5	Cleaning the jet assembly	17
5.6	Replacing fuses	17
5.7	Maintenance Data Sheet	18
6	Technical data	18
7	Display texts and diagrams	19
7.1	Normal mode	19
7.2	The System ready dialog	20
7.3	The syringe selection	20
7.4	The method stack dialog	21
7.5	The parameter type selection dialog	22
7.6	The parameters input dialog	23
7.7	The Setup dialog	24
7.8	Error messages	25
	Declaration of Conformity (DoC)	



#### 1 Introduction

#### 1.1 CAMAG Linomat 5

The CAMAG Linomat 5 is a semiautomatic sample application device for use in qualitative, quantitative and preparative Thin-layer chromatography. The Linomat 5 is shipped with EquiLink, RS232 cable and built in interface for connecting to a computer from which the instrument can be programmed.

The parameters for up to 10 application programs can be stored for stand-alone operation without PC- control.

The Linomat 5 is software controlled, visionCATS/winCATS controls and reports all parameters and steps of the TLC analysis including definition of plate material, sample application, derivatization, development and evaluation. All data are handled in a cGMP/cGLP-compliant manner.

#### 1.2 Precautions


- Please read this operating manual before starting the installation! This manual contains information and warnings the user has to follow to ensure reliable operation of the instrument.
- Some interior parts of the instrument are under AC power. Careless and improper use can cause injury. Unauthorized manipulations can cause damage




- This sign indicates (on instrument and in this manual) that failure to take note of the accompanying information may result in damage of the instrument
- The instrument is manufactured and tested in accordance with the respective European safety publications shown on the Declaration of Conformity (DoC). The instrument complies with safety class 1 and has been designed for indoor use only (IP 20). Further, this device has passed the CAMAG Quality Assurance tests and has been delivered in safe operation condition. For detailed instrument data see chapter technical data
- Attention: For safety reasons the instrument may only be used for the purposes described in the operating manual
- To avoid injury use adequate safety equipment (protective goggles, gloves etc. if applicable) when working with the instrument
- Before first operation, check whether the voltage shown on the instrument matches your local mains voltage. The power cord may only be connected to a grounded, fused (not higher than



- 16A) outlet. Do not use extension cords without ground contact
- When working with the fluids of the instrument, be sure to take the appropriate caution (protect your eyes from direct contact with liquid)
- The instrument may be used only by properly trained laboratory staff
- The instrument may not be used in rooms with danger of explosions
- The instrument contains highly sophisticated electronics and optical parts. It may be operated only in a non-condensing atmosphere in the temperature range outlined in the chapter "Technical Data". Before installation and use, the instrument should be acclimated properly
- Use a damp lint free cloth for cleaning the instrument surface.
   Do not employ aggressive detergents
- Protect yourself and the instrument from electrostatic shock which can cause damage to the electronic parts
- Only authorized personnel may open the instrument. Service and repair is only to be performed by trained specialists. Use spare parts and consumables supplied by CAMAG only. The warranty is voided if parts from other sources are used. Check the service manual before you start service to reduce productspecific risks
- The power cord has to be removed before the instrument is opened. It is not permitted to work on an instrument that has been opened and is connected to the power supply



- Use only the original, with the instrument delivered power cord type
- If the instrument is found to be defective, it must be switched off and steps must be taken to ensure that it cannot be switched on by mistake
- If liquids penetrate the inside of the instrument, the power has
  to be disconnected immediately. Small amounts of liquid can
  be wiped off and/or dried by means of a hairdryer, with larger
  amounts of liquid a service technician has to be called. A test
  of functionality has to be performed in all cases
- Carry out all safety checks and the preventive maintenance as recommended by the manufacturer in order to assure your personal safety and the full functionality of the instrument. Have





- an authorized service specialist perform any service not described by this manual
- See original manufacturers' manuals for further safety data on third party equipment supplied with the system
- Lift/move/transport the system with the necessary care and with sufficient manpower (install the transport security devices if applicable, transport it only in the original packaging)
- The safety of any system incorporate with the equipment is the responsibility of the assembler of the system



This symbol indicates that this equipment must not be disposed
of as unsorted municipal waste but is to be collected separately
as electrical and electronic equipment (WEEE-Directive
2002/96/EC). To properly recycle the instrument or parts of it
you are requested to send the equipment back to the distributor, producer or an adequate collection system at the end of its
life. This will have potential effects on the environment and
human health

### 1.3 Parts supplied

Part no	Units	Description		
695.0014	1	Sample Syringe 100µL		
140.7808-1	1	Tool kit, comprising:		
715.2329	1	Down holder for TLC sheets		
715.2328	1	Down holder for TLC plates		
140.0450-1	1	Pneumatics connecting kit		
125.1021	1	RS232 connecting cable		
362.0006	2	Fuses 1AT		
705.0018	1	Allen wrench 2.5mm		
	1	Respective power cord		
B.022.7808E	1	Instruction manual (English)		

### 1.4 Spare parts

Part no	Description
695.0014	Sample Syringe 100µL
695.0015	Sample Syringe 500µL
022.7830	Plunger for syringe 100µL
022.7831	Plunger for syringe 500µL
115.0570	Venturi jet assembly



### 2 Installation

#### 2.1 Instrument

- Observe the environmental requirements (2.2 Installation environment) when setting up the instrument.
  - ✓ Carefully unpack all components and accessories listed on the shipping list from of the upper part of the packing. Make sure the shipment is complete.

The plate support has been secured with foam pieces for shipment. To unlock the plate support, remove all foam pieces from the slit in the instrument base.




Fig. 1: The instrument

- ✓ Remove the dosage turret carefully from the package and place it on the turret support (fig. 2).
- ✓ Make sure that the dosage turret has been properly positioned on the fixation (3) before tightening the two turret screws (2 and 4) by means of the screwdriver supplied.
- ✓ Tighten screw (1) by hand.
- ► Connect the turret cable with connector (2) as shown in fig. 3.




Fig. 2: Mounting the dosage turret

#### 2.2 Installation environment

The place for installation should meet the following requirements:

Bench space Width 360mm

Depth 510mm (add space for cables approx. 100mm)

Height 410mm Weight 12.5kg

Add adequate space for a PC next to the instrument

Operating The temperature should be within a range of 18 to 35 temperature degrees centigrade and free from significant variations.

Humidity Humidity and temperature conditions must not cause

condensation.

Adequate ventilation free from acidic, alkaline or other Atmospheric conditions

gas that may corrode metal or painted surfaces must

be secured.

#### **Further requirements:**

- Do not place the instrument in a location where the temperature significantly changes (e.g. under an air conditioning duct or by a window). Significant changes in temperature will affect the performance of the unit.
- Do not place the instrument in direct sunlight. Direct sunlight may create significant temperature changes which will affect the performance of the system.
- Do not use the instrument in an environment with moving ambient air (draft).
- Do not expose the instrument to any strong vibration or shock.
- Avoid placing the instrument near equipment that radiates heat. Do not place the instrument near gas burners, electric heaters or ovens.



- Do not place the instrument near equipment that generates intense magnetic fields such as electric welding equipment, high frequency furnaces, pole transformers, etc.
- Protect the instrument from excessive dust.
- Connect the instrument to power lines that are free from sudden changes or voltage fluctuations.
- If you must use power motor driven equipment (such as a stirrer or shaker) in the same line as your instrument, ensure that a noise reduction unit is in the same power line

#### 2.3 Conditions for the installation

Confirm that the following requirements exist before installing the instrument:

Power supply and ground

Line voltage: 100 - 240V (see rating plate on instrument).

Frequency: 50 / 60 Hz Power capacity: 20 W

Ground A grounded outlet should be within 2 meters of

terminal: the instrument.

Gas supply

N2 or Required pressure 4 - 6 bar (58 -87 psi)

compressed air Gas consumption 1 L/min

### 2.4 Gas supply

- ✓ Connect the supplied pressure tubing to the gas inlet nipple (fig. 3, item 4) on the backside of the instrument.
- ✓ Connect the other end of the tubing to the gas supply.

If the outlet of your gas supply does not fit, the  $N_2$ /compressed air connection kit allows connection to standard  $\frac{1}{4}$ " tubing.

- ✓ Adjust the primary pressure to 4 6 bar (58-87 psi).
- The message: "WAITING FOR SYSTEM PRESSURE" will be displayed if the instrument is started with too low pressure.
- ► If long thin tubing is used it is recommended to increase the primary pressure to 6-8 bar (87-115 psi).

#### 2.5 Mains connection

- ✓ Ensure that the voltage shown on the rating plate (left side) matches that of the mains.
- ✓ Use the supplied power cord for connecting the instrument to the mains outlet.




Fig. 3: Connections. 1 = RS232 Interface, 2 = dosage turret connection, 3 = Power switch and mains connection with fuse box, 4 = gas connection

#### 2.6 Connection to the PC

The CAMAG Linomat 5 is shipped with a 9-pin interface cable for connection to the PC.

✓ Connect the cable to the instrument and to a free COM-port of the PC to enable data transfer.

### 3 Getting started

### 3.1 Automatic instrument monitoring

A sensor in the CAMAG Linomat 5 monitors the pressure of the gas supply and generates an error message if the pressure drops below 3.5 bar during spray-on application.

All valves are controlled in such way that they close if the power fails or the instrument is switched off.

### 3.2 Checking system functions

✓ Switch on the instrument.

The Linomat 5 will now perform the following checks and tasks:

- 1. Initializing the plate stage drive
- 2. Initializing the syringe drive.

If startup was successful the display will read:

\*\* Linomat 5 \*\*
Instrument ready

In this phase, only physical blocking of a drive can cause an error. If necessary, check if all transport locks have been removed and whether there are foreign objects in the instrument. (e.g. below the plate stage).



### 3.3 Inserting a dosage syringe

#### Installation of a syringe (see fig. 4):

- ✓ Make sure plate stage is in front position and the syringe lever arm is at the top of the turret.
- ✓ Fill the syringe with at least the minimum amount indicated on the display.
- ✓ Use the left hand to move the lever arm (1) to the left and hold it there.
- ✓ Carefully insert the syringe through the syringe holder (2) into the spray head (3) and lower the syringe carefully until it rests in the spray head. This way the needle is guided through the spray nozzle.
- The volume of the syringe in use must be set in the user dialog. If a wrong syringe volume is entered, the application volumes will all be wrong.




Fig. 4: Installation of the syringe. 1 = Lever arm, 2 = syringe holder, 3 = Spray head assembly with nozzle, 4 = slide coupler, 5 = scale pointer

#### **Entering syringe volume:**

- ✓ Press the **DIALOG** key on the instrument to start the user dialog.
- ✓ Use the arrow keys and to select the volume of the syringe in use.
- ✓ Confirm by pressing the **ENTER** key on the instrument.



### 3.4 Setting the application position

The application position Y is set by moving the dosage turret to the proper position.

- ✓ Press the slide coupler (4) backwards till the lock pin is released.
- ✓ Move the turret to the desired application position and then release the slide coupler. You can read the position (5) on the ruler at the front edge of the turret bridge.

### 3.5 Loading a TLC plate or foil

All plate sizes held by both front and rear ledge (20x20, 20x10):

- ✓ Move the plate lift lever (1) situated on the right side of the plate stage to the rear position in order to lower the table.
- ✓ Now flip the front ledge (2) of the stage down (forward) making the plate table freely accessible.
- ✓ Place the plate (3) in the rear left corner of the table so that it touches the rear stage ledge and both positioning pins (4).
- ✓ Now flip up the front ledge (2) and move the plate lift lever (1) carefully to the forward position. Guide the lever to avoid that the stage moves up too quickly.




Fig. 5: Plate stage. 1 = plate stage lever, 2 = front ledge, 3 = TLC-plate, 4 = pin.

All sizes not held by both front and rear ledge (10x10 or smaller):

- ✓ Position the plate to be applied on as described above
- ✓ Position another plate of same thickness under the front ledge for thickness compensation.
- ✓ Position the plate holder behind the TLC plate to add stability.




Fig. 6: Plate stage for (2x)10 x 10 cm plates.

#### TLC foils (1):

Aluminium foils may not lay flat on the stage but can be held in place by means of the down holder for foils (2). It is employed according to fig. 7 and may be left in place if only foils are being used.

The positioning of foils is similar to that of TLC-plates although the edge of the foil must be placed under the down holder and in contact with the positioning pins.




Fig. 7: Plate stage with foil. 1 = foil, 2 = down holder for foils.

### 3.6 Waste plate

To ensure proper spraying all newly filled syringes are started by predosing a small volume at the waste position. The waste position is below the spray nozzle wherever the turret is currently positioned.

- ✓ Put a 20x10cm TLC plate (waste plate) on the instrument housing below the spray nozzle.
- To avoid predosing into the slit in the housing make sure not to position the turret between (approx.) 40 and 65mm (3).
- Move the waste plate from time to time so the spray does not always hit the same position.




Fig. 8: 1 and 2 = possible positions for waste plate, 3 = no spray area

### Manual operation

### 4.1 Display and function of keys

#### **Control lamps**

The control lamps in the display indicate the following status:

- POWER ON: this lamp is lit when the instrument is switched on.
- ON LINE: this lamp is lit when the instrument is ON-LINE with winCATS or is communicating with visionCATS

#### The following keys are available:

#### **DIALOG**

Selects the syringe size selection display (System ready) or activates the manual parameter input mode (Method stack).

#### **ARROW KEYS ↑** and **▼**

Select the next / previous parameter (or method) in the dialog.

#### **ENTER**

Confirms the current parameter, switches to the selected dialog or (re-)starts the (last) selected method.

#### **END**

Interrupts the current application or switches back to the previous dialog.

#### **RESET**

Resets and re-initializes the instrument. The current application is interrupted but all parameters remain in memory.

#### **RUN**

Shows the selection of stored methods.

✓ Use the arrow keys and to select a method to start.



✓ Press **ENTER** to start the selected method or **DIALOG** to change the parameters of the method locally.

### 4.2 Stand-alone operation

The instrument can be used in stand-alone mode. Up to 10 methods can be stored and used locally.

Note that locally started methods cannot be used in an analysis in visionCATS/winCATS and that local parameter changes are not documented in the software.

#### **Procedure for OFF-LINE start:**

- ✓ Prepare the instrument for sample application.
- ✓ Press the RUN key
   The most recently used method is displayed. With the and
   ✓ keys you can select another locally stored method.
- ✓ Press the ENTER key to start the displayed method. The last method is re-started if ENTER is pressed in the "System ready" mode.

#### Procedure for local input of parameters:

- ✓ Press the RUN key The most recently used method is displayed. With the or vert key you can select another locally stored method.
- ✓ Press the **DIALOG** key to locally change parameters of the selected method.
  - Press the ♠ or ▼ key to select a parameter group (global parameters, layout parameters, track assignment parameters, save method) for local input and press **ENTER** to confirm.
- ✓ Now you can step through the parameters of the group by means of the ENTER key or use the ^ or ▼ key to change the current parameter. Press ENTER to confirm the new value and move on the next parameter.

#### Procedure for saving a changed method:

- ✓ Select the parameter group Save method or press the END key after having changed the parameters.
- ✓ Select the memory position to save the new method in by means of the or key. Press ENTER to confirm position and execute the save routine.

#### Procedure for local delete of a saved method:

✓ Press the RUN key
The most recently used method is displayed. With the or vert key you can select the locally stored method to delete.



- ✓ Press the key combination RUN-END (press and hold RUN, then press END) and confirm deleting with ENTER.
- Flow diagrams for these dialogs are depicted in part 7 Display texts.

### 4.3 Aborting an application in progress

✓ Press the END key The application in progress is aborted. An aborted application cannot be continued.

### 4.4 Resetting the instrument

#### **Procedure:**

✓ Press the RESET key This will abort all activities in progress and re-initialize the instrument.

### 4.5 Docking the syringe

After a method has been started the display will show the required volume and the vial or ID of the sample to be used.

- ✓ Draw up the required volume of the sample indicated.
- ✓ Mount the syringe as described in chapter 3.3
- ✓ Press ENTER

The syringe drive will now perform an automatic docking at the indicated volume, then the spraying is started at the waste position and finally sample application starts.

If you do not wish to perform automatic docking after mounting the syringe, you can move the syringe drive manually by means of the ▼ key. Press both arrow keys to move fast.



### Maintenance and trouble shooting

#### 5.1 General maintenance

The instrument may be serviced only by authorized technicians familiar with its technical/functional properties. In daily use the Linomat 5 requires maintenance about once a year. The chapter Maintenance Data Sheet informs about the parts and their frequency to be changed.

### 5.2 The Setup dialog

By means of the Setup dialog you can select the communication baud rate between the instrument and the PC, the language and LCD contrast.

#### To enter the Setup dialog:

- ✓ Hold down the **DIALOG** key and press the **RESET** key.
- ✓ Release the **RESET** key, then release the **DIALOG** key.

The display will now show **LINO5 Vx.xx.xx Setup Mode** where x.xx.xx is the actual firmware version of the Linomat 5.

✓ Press the **ENTER** key to start the Setup dialog.

#### Set baud rate

Start the Setup dialog according to the procedure above. The display now shows the current baud rate (default is 19'200).

- ✓ Press the or vert key to change the setting.
- ✓ Press the ENTER key to accept the current selection (and move on to the language selection).
- ✓ Press the **RESET** key to quit the user dialog.

### Set display language

Start the user dialog according to the procedure above. The display now shows the current baud rate.

✓ Press the ENTER key to move on to the language selection display.

The current language setting is displayed.

- ✓ Press the  $\blacktriangle$  or  $\blacktriangledown$  key to change the setting.
- ✓ Press the ENTER key to accept the current selection (and move on to the LCD contrast selection).
- ✓ Press the **RESET** key to quit the user dialog.



#### Set the LCD contrast

Start the user dialog according to the procedure above. The display now shows the current baud rate.

✓ Press the ENTER key twice to move on to the LCD contrast selection display.

The current contrast setting is displayed.

- ✓ Press the  $\blacktriangle$  or  $\blacktriangledown$  key to change the setting.
- ✓ Press the **ENTER** key to accept the current selection.
- ✓ Press the **RESET** key to restart the instrument.

### 5.3 Handling of the HAMILTON-syringe

To achieve a maximum lifetime for the gas-tight HAMILTONsyringe the following points must be observed:

- Ensure extremely clean conditions when handling the syringe body and plunger.
- Never lubricate the plunger.
- The Teflon tip of the plunger is very sensitive. Even slightest scratches can cause a leak.
- Plungers and syringe bodies are exchangeable.
- To fill the sample dosage syringe make sure the needle is not pushed against the bottom of vial and that no air bubbles remain in the needle.
- The plunger must not be exposed to boiling water.

The following replacement plunger can be obtained:

022.7830 Replacement plunger for 100  $\mu$ L syringe

022.7831 Replacement plunger for 500  $\mu L$  syringe

### 5.4 Adjusting the spray gas

Spray gas flow is controlled by a gas pressure regulator, factory adjusted to 1.0 bar. It may be suitable to increase this pressure to achieve better spray quality when using highly viscous liquids.

#### **Procedure:**

✓ Turn the spray pressure regulator (1) clockwise to increase pressure. The pressure can be read from an external gauge you can attach at (2). Adjust pressure to 1.5 bar.




Fig. 9: 1 = pressure regulator, 2 = connection for external gauge, 3 = gas supply,

Change the pressure <u>only while sample application is in progress</u> to get a correct reading from the gauge.

### 5.5 Cleaning the jet assembly

The jet assembly, a part of spray nozzle, consists of a square outlet into which a circular needle is inserted. At the four free corners gas flows at high speed and sprays the sample emerging from the needle tip onto the TLC plate. Clogged corners make the jet stream irregular.

✓ Unscrew the jet assembly and wash it with a suitable solvent (ultrasonic bath).

### 5.6 Replacing fuses

#### **Procedure:**

- ✓ First disconnect the power cord from the mains connection of the instrument.
- ✓ Pull out the fuse holder next to the power connection with a small screwdriver.
- ✓ Replace the fuse(s). The label 250 V 1.0 AT (slow) is valid to the full voltage range.
- ✓ Push back the fuse box.
- ✓ Reconnect the power cord.




Fig. 10: Mains connection, switch and fuse holder



### 5.7 Maintenance Data Sheet

# CAMAG Maintenance Data Sheet

### Linomat 5

February 2015

Purpose	The maintenance data sheet informs about maintenance interval of the respective
	instrument as well as the proposal for IQ/OQ interval if applicable. In addition, it
	identifies consumable parts with the respective replacement cycle.

Maintenance interval			
Maintenance	12 Months		
IQ/OQ	12 Months		

Consumable parts			
Part No.	Description	Replacement cycle	
695.0014	Syringe 100µl	24 Months	
695.0015	Syringe 500µl	24 Months	
115.0570	Spray nozzle	24 Months	
	Plunger for 100µl syringe	3 Months	
022.7831	Plunger for 500µl syringe	3 Months	

### 6 Technical data

Power connections	100 – 240 V 50 / 60 Hz 20 W			
Dosing syringe	Choice of 100 μL or 500 μL, gas-tight			
Dosing volume	0,1μL to 5000μL (i.e. 50 fillings of 100μL syringe per application!)			
Syringe drive	Stepper motor 1600 steps / rev.:			
	100nL = 120 steps for 100μL syringe			
	100nL = 24 steps for 500μL syringe			
Stage drive	Stepper motor 3200 steps/rev. 8 steps = 0.1 mm			
	Maximum speed 250 mm/s with acceleration ramp			
	Positions programmable: 5.0 – 195.0 mm in 0.1 mm increment			
In- output connec-	RS232			
tion				
Measures	Width 360mm, depth 510mm, height 410mm			
Weight 12.5 kg				




## Display texts and diagrams

### 7.1 Normal mode

Text	Description		
FW V: 1.00.01 *****	Firmware version at start up Device Power on		
** CAMAG ** **LINOMAT5**	The system was properly initialized and is waiting for an application.		
LINOMAT5 SYSTEM READY	The system was properly in initialized and is waiting for application.		
THIS METHOD? TEST 87364	Selection of a stored method		
WAITING FOR SYSTEM PRESSURE	System pressure too low (46bar). Device is waiting for gas.		
GO TO PARK POSITION	Plate table and syringe lever are moved to the parking position.		
INIT TABLE DRIVE	Initialization of specified drive.		
<b> </b>	Predosage is running for the 2. of 12 samples and the 5. of total 32 bands		
<b> </b>	Dosage is running for the 3. of total 5 samples and the 9. of total 18 bands. 4.7µl of the 2. of total 4 bands from sample 3 has been applied.		
OPERATION COMPLETED	Application has been properly completed		
LINOMAT 5 TESTING DRIVES	Drives test is running		
** PAUSE! ** 12 min 48 s	Instrument is in pause mode since 12 min and 48 sec.		
** ERROR 38 ** TEXT TEXT TEXT	See chapter 7.4 for error messages		



### 7.2 The System ready dialog



### 7.3 The syringe selection


To enter the Syringe selection, press the **DIALOG** key when the instrument display shows "System ready".

Text	Description
SELECT SYRINGE	Select syringe volume with UP and DOWN
100/500ul	keys



### 7.4 The method stack dialog


To enter the Method selection dialog (Method stack), press the **RUN** key when the instrument display shows System ready.





### 7.5 The parameter type selection dialog


To enter this parameter type selection dialog, press the **DIALOG** key when the instrument display shows the name of the method to be edited (in the method stack dialog).





### 7.6 The parameters input dialog


To enter this parameter input dialog, press the **ENTER** key when the display shows the name of the parameter type to be edited (in the parameter type selection dialog).





### 7.7 The Setup dialog

To enter the Setup dialog, hold down the dialog key wile switching on the instrument.





### 7.8 Error messages

Error#	Description		
30	Pneumatic pressure too low at startup		
31	Wrong syringe mounted.		
	The syringe specified in the method does not match the syringe currently registered in the Linomat 5. To register another syringe: press the <b>DIALOG</b> key ion the System ready dialog.		
33	Combination errors in method parameters during local		
	parameter input, see also error 36.		
34	Dosage volume 0µl not possible		
35	Dosage volume >100/500ul not allowed for one band. The Linomat 5 cannot split the volume of one band on several syringe fillings.		
36	Selected dosage parameters unrealistic:		
	The overspray factor (OF) is calculated as:		
	Dosage time=Application volume/Dosage speed		
	Band time=Band length/band velocity		
	OF=Dosage time/Band time		
	The band velocity is automatically optimized between 1		
	and 20mm/s (target is 10mm/s).		
	Possible errors:		
	1 The overspray factor (OF) is close to zero.		
	2 The max. distance sprayed is > 50m (OFxBand length)		
37	Selected dosage parameters unrealistic:		
20	Turret drive speed too high!		
38	Selected dosage parameters are unrealistic:		
55	Turret drive speed too low!  Method storage empty!		
56	Content of activated method incorrect (data loss )		
59			
100	Interrupt of method by main interrupt Syringe drive: end switch (LSX) active		
100	Stage drive: end switch (LSX) active		
104			
104	Syringe drive: end switch was expected Stage drive: end switch was expected:		
105	Syringe drive: initialization has not been executed		
107	Stage drive: initialization has not been executed		
120	Pressure loss during dosage		
	5 5		
121	Internal device voltage 24V got lost or too low		

CAMAG (Switzerland) · Sonnenmattstrasse 11 · CH-4132 Muttenz 1 Telephone +41 61 467 34 34 · Fax +41 61 461 07 02 · E-Mail: <u>info@camag.com</u>

CAMAG (Germany) · Bismarckstraße 27-29 · D-12169 Berlin Telephone +49 30 516 55 50 · Fax +49 30 795 70 73 · E-Mail: infoberlin@camag.com

CAMAG Scientific Inc. (USA) · 515 Cornelius Harnett Drive · Wilmington, NC 28401 Telephone 800 334 3909 · Fax 910 343 1834 · E-Mail: tlc@camag.com

www.camag.com



### **EC – Declaration of Conformity**

We, CAMAG Chemie-Erzeugnisse und Adsorptionstechnik AG Sonnenmattstrasse 11 4132 Muttenz Switzerland

declare under our sole responsibility that the product

#### **CAMAG® Linomat 5**

Product name

022.7808

Article number(s)

to which this declaration relates is in conformity with the following provisions of directive(s):

- 2006/95/EC
- 2004/108/EC

Following standard(s) or other normative document(s):

EN61010-1: 2010EN61326-1: 2013

Muttenz, 26 March 2015

Johek-Ralme

Walter Rahm, Head of Quality Management

SWISS HADE

CE



Notes:			
		 	<u></u> .
		 	<del></del>
,		 	
<del></del>		 	
	<del></del>		

Notes:			
		 	<u></u> .
		 	<del></del>
,		 	
<del></del>		 	
	<del></del>		