

## Blast plasma freezer

### Applications

- Designed for rapid freezing of human plasma



### Spécifications

- Maintaining temperature at -40°C at the end of the freezing cycle
- Navigation by touch screen
- Security HP pressostat
- Performances :
  - 11°C to the core within 60 minutes
  - 30°C to the core within 90 minutes

### Avantages

- Setting temperature parameters protected by password
- Launch cycle by push-button
- Display and digital recording of the temperature
- Data storage in a non-volatile memory flash
- Data logging
- Control freezing functioning by a temperature differential T° / Temps.



Reference	Name	280 ml blood bag capacity	Condenser		External	
			Air	Water	Condenser	Refrigeration Unit
CRP20A	CRP 20 / A	20	X			
CRP40A	CRP 40 / A	40	X			
CRP40E	CRP 40 / E	40		X		
CRP40AGD	CRP 40 / A / GD	40	X			X
CRP40ACD	CRP 40 / A / CD	40	X		X	
CRP100E	CRP 100 / E	100		X		X
CRP100AGD	CRP 100 / A / GD	100	X			X

### ACCESSORIES

Reference	Models	Description
CRP/CHAR	CRP 40 CRP 100	Loading trolley on castors with racks and shelves
CRP/RACK	CRP 40 CRP 100	Racks and shelves

### OPTIONS

Reference	Models	Description
CRP:INSTAL	CRP 40 CRP100	Installation on demand